

CYP11B1 Polyclonal Antibody



Catalog Number: 162415

Product name

CYP11B1 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to CYP11B1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 234-503 of human CYP11B1 (NP_000488.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:200

Validations



Immunofluorescence - CYP11B1 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using CYP11B1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.