

# SLC17A6 Polyclonal Antibody



Catalog Number: 162493

## Product name

SLC17A6 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to SLC17A6

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 500 to the C-terminus of human SLC17A6 (NP\_065079.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

## Validations



Western blot - SLC17A6 Polyclonal Antibody

Western blot analysis of extracts of SH-SY5Y cells, using SLC17A6 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.