

CSNK1D Polyclonal Antibody



Catalog Number: 162753

Product name

CSNK1D Polyclonal Antibody

Specificity

Rat

Antibody description

Polyclonal antibody to CSNK1D

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 250-350 of human CSNK1D (NP_620693.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - CSNK1D Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using CSNK1D antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.