

BMP1 Polyclonal Antibody



Catalog Number: 163158

Product name

BMP1 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to BMP1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 725-986 of human BMP1 (NP_006120.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:200

Validations



Immunofluorescence - BMP1 Polyclonal Antibody
Immunofluorescence analysis of HeLa cells using BMP1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.