

NEFM Polyclonal Antibody

Catalog Number: 163181

Product name

NEFM Polyclonal Antibody

Specificity

Mouse, Rat

Antibody description

Polyclonal antibody to NEFM

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 200-300 of human NEFM (NP_005373.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

Validations




Western blot - NEFM Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using NEFM antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunofluorescence - NEFM Polyclonal Antibody

Immunofluorescence analysis of mouse brain using NEFM antibody at dilution of 1:100. Blue:



DAPI for nuclear staining.

NEFM Polyclonal Antibody



Catalog Number: 163181

Immunofluorescence - NEFM Polyclonal Antibody

NEFM antibody at dilution of 1:100. Blue: DAPI for nuclear staining.

Immunofluorescence analysis of rat brain using