

CD70 Polyclonal Antibody



Catalog Number: 163350

Product name

CD70 Polyclonal Antibody

Specificity

Mouse,Rat

Antibody description

Polyclonal antibody to CD70

Preparation

Antigen: A synthetic peptide of human CD70.

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

IF

Dilutions

IF 1:50 - 1:200

Validations

Immunofluorescence - CD70 Polyclonal Antibody

Immunofluorescence analysis of rat spleen using CD70 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.