

# NPAS2 Polyclonal Antibody



Catalog Number: 163377

## Product name

NPAS2 Polyclonal Antibody

## Specificity

Human, Rat

## Antibody description

Polyclonal antibody to NPAS2

## Preparation

Antigen: A synthetic peptide of human NPAS2.

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:100

## Validations



Western blot - NPAS2 Polyclonal Antibody

Western blot analysis of extracts of A-431 cells, using NPAS2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.