

CCL22 Polyclonal Antibody



Catalog Number: 163638

Product name

CCL22 Polyclonal Antibody

Specificity

Human, Rat

Antibody description

Polyclonal antibody to CCL22

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-93 of human CCL22 (NP_002981.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions


IHC 1:50 - 1:100

Validations



Immunohistochemistry - CCL22 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat kidney using CCL22 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - CCL22 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon using CCL22 antibody at dilution of 1:100 (40x lens).