

APOBEC3F Polyclonal Antibody



Catalog Number: 163810

Product name

APOBEC3F Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to APOBEC3F

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-65 of human APOBEC3F (NP_660341.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Immunohistochemistry - APOBEC3F Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat ovary using APOBEC3F antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - APOBEC3F Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat heart using APOBEC3F antibody at dilution of



1:100 (40x lens).


Immunohistochemistry - APOBEC3F Polyclonal Antibody

APOBEC3F Polyclonal Antibody



Catalog Number: 163810


Immunohistochemistry of paraffin-embedded human lung using APOBEC3F antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - APOBEC3F Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human lung cancer using APOBEC3F antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - APOBEC3F Polyclonal Antibody


Immunohistochemistry of paraffin-embedded

human stomach using APOBEC3F antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - APOBEC3F Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using APOBEC3F antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - APOBEC3F Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using APOBEC3F antibody at dilution of 1:100 (40x lens).