

# CCL25 Polyclonal Antibody



Catalog Number: 163856

---

## Product name

CCL25 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to CCL25

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 24-150 of human CCL25 (NP\_005615.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IF

## Dilutions

IF 1:50 - 1:200

## Validations



Immunofluorescence - CCL25 Polyclonal Antibody

Immunofluorescence analysis of MCF7 cells using CCL25 antibody . Blue: DAPI for nuclear staining.