

# TP73 Polyclonal Antibody



Catalog Number: 163893

## Product name

TP73 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to TP73

## Preparation

Antigen: Recombinant protein of human TP73

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:200 - 1:1000

## Validations



## Western blot - TP73 Polyclonal Antibody

Western blot analysis of extracts of MDA-MB231 cells, using TP73 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.