

CCL8 Polyclonal Antibody



Catalog Number: 163976

Product name

CCL8 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to CCL8

Preparation

Antigen: A synthetic peptide of human CCL8

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB

Dilutions

WB 1:200 - 1:1000

Validations

Jurkat



Western blot - CCL8 Polyclonal Antibody

Western blot analysis of extracts of Jurkat cells, using CCL8 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.