

PACSIN2 Polyclonal Antibody



Catalog Number: 164060

Product name

PACSIN2 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to PACSIN2

Preparation

Antigen: A synthetic peptide of human PACSIN2

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:200 - 1:1000

IHC 1:50 - 1:200

IF 1:20 - 1:50

Validations



Western blot - PACSIN2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using PACSIN2 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.