

# PRIM2 Polyclonal Antibody



Catalog Number: 164260

## Product name

PRIM2 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to PRIM2

## Preparation

Antigen: A synthetic Peptide of human PRIM2

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:1000

## Validations



Western blot - PRIM2 Polyclonal Antibody

Western blot analysis of extracts of NCI-H460 cells, using PRIM2 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.