

ABAT Polyclonal Antibody



Catalog Number: 164471

Product name

ABAT Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to ABAT

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 29-300 of human ABAT (NP_001120920.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations




Western blot - ABAT Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ABAT antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.



Immunohistochemistry - ABAT Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat kidney using ABAT antibody at dilution of 1:100




(40x lens).

ABAT Polyclonal Antibody

Catalog Number: 164471

Immunohistochemistry - ABAT Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human liver injury using ABAT antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - ABAT Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human gastric cancer using ABAT antibody at



dilution of 1:100 (40x lens).

Immunofluorescence - ABAT Polyclonal Antibody

Immunofluorescence analysis of A549 cells using ABAT antibody . Blue: DAPI for nuclear staining.