

POLD1 Polyclonal Antibody



Catalog Number: 164477

Product name

POLD1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to POLD1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 808-1107 of human POLD1 (NP_002682.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations



Western blot - POLD1 Polyclonal Antibody


Western blot analysis of extracts of various cell lines, using POLD1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Basic Kit .Exposure time: 90s.

Immunofluorescence - POLD1 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using POLD1 antibody . Blue: DAPI for nuclear staining.




Immunofluorescence - POLD1 Polyclonal Antibody

POLD1 Polyclonal Antibody

Catalog Number: 164477

Immunofluorescence analysis of A549 cells using



POLD1 antibody.

Immunofluorescence - POLD1 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using POLD1 antibody . Greenuff1agamma H2AX (phospho S139) used for DNA damage marker. Blue: DAPI for nuclear staining. Cells were treated by UV-A laser.