

TYROBP Polyclonal Antibody



Catalog Number: 164501

Product name

TYROBP Polyclonal Antibody

Specificity

Mouse, Rat

Antibody description

Polyclonal antibody to TYROBP

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-113 of human TYROBP (NP_003323.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions


IHC 1:50 - 1:200

Validations



Immunohistochemistry - TYROBP Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat heart using TYROBP antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - TYROBP Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using TYROBP antibody at dilution of 1:100 (40x lens).