

ZEB2 Polyclonal Antibody

Catalog Number: 164601

Product name

ZEB2 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to ZEB2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1050-1214 of human ZEB2 (NP_055610.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations




Western blot - ZEB2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ZEB2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - ZEB2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using ZEB2 antibody at dilution of 1:100




(40x lens).

ZEB2 Polyclonal Antibody

Catalog Number: 164601

Immunohistochemistry - ZEB2 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human breast using ZEB2 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - ZEB2 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded mouse liver using ZEB2 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - ZEB2 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded mouse testis using ZEB2 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - ZEB2 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded mouse spinal cord using ZEB2 antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - ZEB2 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded mouse brain using ZEB2 antibody at dilution of



1:100 (40x lens).

Immunofluorescence - ZEB2 Polyclonal Antibody

Immunofluorescence analysis of C6 cells using ZEB2 antibody at dilution of 1:100. Blue: DAPI for



nuclear staining.

Immunofluorescence - ZEB2 Polyclonal Antibody

Immunofluorescence analysis of NIH/3T3 cells

ZEB2 Polyclonal Antibody



Catalog Number: 164601

using ZEB2 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.