

GATA3 Monoclonal Antibody



Catalog Number: 164605

Product name

GATA3 Monoclonal Antibody

Specificity

Human

Antibody description

Monoclonal antibody to GATA3

Preparation

Antigen: A synthetic peptide of human GATA3

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG

Applications

WB, IF, FC


Dilutions

WB 1:500 - 1:2000

IF 1:20 - 1:50

FC 1:20 - 1:50

Validations



Western blot - GATA3 Monoclonal Antibody

Western blot analysis of extracts of various cell lines, using GATA3 antibody. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.