

# LIPA Polyclonal Antibody



Catalog Number: 164846

## Product name

LIPA Polyclonal Antibody

## Specificity

Human, Rat

## Antibody description

Polyclonal antibody to LIPA

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 110-399 of human LIPA (NP\_001121077.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200

## Validations



Western blot - LIPA Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using LIPA antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.