

RBP2 Polyclonal Antibody

Catalog Number: 164996

Product name

RBP2 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to RBP2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-134 of human RBP2 (NP_004155.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:100

Validations




Western blot - RBP2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using RBP2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.



Immunohistochemistry - RBP2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver injury using RBP2 antibody at dilution




of 1:100 (40x lens).

RBP2 Polyclonal Antibody

Catalog Number: 164996

Immunohistochemistry - RBP2 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human stomach using RBP2 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - RBP2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using RBP2 antibody at dilution of 1:100



(40x lens).

Immunofluorescence - RBP2 Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using RBP2 antibody . Blue: DAPI for nuclear staining.