

CCL8 Polyclonal Antibody



Catalog Number: 165122

Product name

CCL8 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to CCL8

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 23-99 of human CCL8 (NP_005614.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100


IF 1:50 - 1:100

Validations



Immunohistochemistry - CCL8 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human liver using CCL8 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - CCL8 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human prostate cancer using CCL8 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - CCL8 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using CCL8 antibody at dilution of

CCL8 Polyclonal Antibody



Catalog Number: 165122

1:100 (40x lens).