

DPP3 Polyclonal Antibody



Catalog Number: 165163

Product name

DPP3 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to DPP3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human DPP3 (NP_569710.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations




Western blot - DPP3 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using DPP3 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - DPP3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat spinal cord using DPP3 Antibody at dilution of



1:100 (40x lens).

DPP3 Polyclonal Antibody



Catalog Number: 165163

Immunohistochemistry - DPP3 Polyclonal Antibody spleen using DPP3 Antibody at dilution of 1:100
Immunohistochemistry of paraffin-embedded rat (40x lens).