

GRIN1 Polyclonal Antibody

Catalog Number: 165228

Product name

GRIN1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to GRIN1

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 600-700 of human GRIN1 (NP_015566.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

Validations




Western blot - GRIN1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using GRIN1 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence - GRIN1 Polyclonal Antibody

Immunofluorescence analysis of rat brain using GRIN1 antibody at dilution of 1:100. Blue: DAPI



for nuclear staining.

Immunofluorescence - GRIN1 Polyclonal Antibody

GRIN1 Polyclonal Antibody



Catalog Number: 165228

Immunofluorescence analysis of SKOV3 cells

using GRIN1 antibody at dilution of 1:100. Blue:
DAPI for nuclear staining.