

OPRM1 Polyclonal Antibody



Catalog Number: 165281

Product name

OPRM1 Polyclonal Antibody

Specificity

Mouse, Rat

Antibody description

Polyclonal antibody to OPRM1

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 300 to the C-terminus of human OPRM1 (NP_000905.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:1000

Validations



Western blot - OPRM1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using OPRM1 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.