

EGR1 Polyclonal Antibody



Catalog Number: 165282

Product name

EGR1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to EGR1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human EGR1 (NP_001955.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC

Dilutions


WB 1:200 - 1:1000

IHC 1:20 - 1:200

IF 1:10 - 1:100

IP 1:20 - 1:50

Validations



Western blot - EGR1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using EGR1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.