

# ADCY3 Polyclonal Antibody



Catalog Number: 165579

## Product name

ADCY3 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to ADCY3

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human ADCY3 (NP\_004027.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:1000

## Validations



## Western blot - ADCY3 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ADCY3 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 60s.