

EYA3 Polyclonal Antibody



Catalog Number: 165595

Product name

EYA3 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to EYA3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human EYA3 (NP_001981.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations




Western blot - EYA3 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using EYA3 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.



Immunohistochemistry - EYA3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat kidney using EYA3 antibody at dilution of 1:100



(40x lens).


EYA3 Polyclonal Antibody



Catalog Number: 165595

Immunohistochemistry - EYA3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human kidney using EYA3 antibody at dilution of




1:100 (40x lens).

Immunohistochemistry - EYA3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded

human gastric using EYA3 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - EYA3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse liver using EYA3 antibody at dilution of 1:100 (40x lens).