

NEDD4L Polyclonal Antibody

Catalog Number: 165697

Product name

NEDD4L Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to NEDD4L

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 836-955 of human NEDD4L (NP_056092.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations




Western blot - NEDD4L Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using NEDD4L antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - NEDD4L Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat kidney using NEDD4L antibody at dilution of



1:200 (40x lens).


NEDD4L Polyclonal Antibody



Catalog Number: 165697

Immunohistochemistry - NEDD4L Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human prostate using NEDD4L antibody at



dilution of 1:200 (40x lens).

Immunofluorescence - NEDD4L Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using NEDD4L antibody at dilution of 1:100. Blue: DAPI for nuclear staining.