

TIMM10B Polyclonal Antibody



Catalog Number: 165913

Product name

Polyclonal

TIMM10B Polyclonal Antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse

Applications

Antibody description

WB, IHC

Polyclonal antibody to TIMM10B

Dilutions

Preparation

WB 1:500 - 1:2000

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-103 of human TIMM10B (NP_036324.1).

IHC 1:50 - 1:100

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Validations

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Western blot - TIMM10B Polyclonal Antibody

