

USP21 Polyclonal Antibody



Catalog Number: 166084

Product name

USP21 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to USP21

Preparation

Antigen: A synthetic peptide of human USP21

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:200 - 1:2000

IHC 1:20 - 1:200

Validations



Western blot - USP21 Polyclonal Antibody

Western blot analysis of extracts of RAW264.7 cells, using USP21 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.