

# NUP98 Monoclonal Antibody



Catalog Number: 166231

## Product name

NUP98 Monoclonal Antibody

## Specificity

Human

## Antibody description

Monoclonal antibody to NUP98

## Preparation

Antigen: Recombinant protein of human NUP98

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG

## Applications


WB, IP

## Dilutions

WB 1:500 - 1:2000

IP 1:20 - 1:50

## Validations



Western blot - NUP98 Monoclonal Antibody

Western blot analysis of extracts of various cell lines, using NUP98 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 10s.