

# Phospho-ULK1-S757 pAb



Catalog Number: 166508

## Product name

Phospho-ULK1-S757 pAb

## Specificity

Human

## Antibody description

Polyclonal antibody to Phospho-ULK1-S757

## Preparation

Antigen: A synthetic phosphorylated peptide around S757 of human ULK1 (NP\_003556.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



### Western blot - Phospho-ULK1-S757 pAb

Western blot analysis of extracts of A431 cells, using Phospho-ULK1-S757 pAb at 1:1000 dilution. A431 cells were treated by EGF (100 ng/ml) at 37°C for 30 minutes after serum-starvation overnight. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% BSA. Detection: ECL Basic Kit. Exposure time: 10s.