

Anti-Spike-RBD Human mAb



Catalog Number: 166656

Product name

Anti-Spike-RBD Human mAb

Specificity

Human

Antibody description

human monoclonal antibody against Spike protein receptor binding domain

Formulation

PBS(pH7.2)

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Human IgG


Applications

ELISA

Dilutions

Paired antibodies to detect SARS-COV-2 S-RBD using the double-antibody sandwich method with antibody(Cat#: 166655)

Validations



Double antibody sandwich Elisa using two anti-S1-RBD antibodies and S1 protein

The two anti-S1-RBD human monoclonal antibodies(Cat#: 166655 and 166656) can simultaneously and stably bind to different epitopes of the SARS-COV-2 S1-RBD protein, with an EC50 value of 0.23 nM.