

# Anti-N protein Human mAb



Catalog Number: 166657

## Product name

Anti-N protein Human mAb

## Specificity

Human

## Antibody description

human monoclonal antibody (IgG) against SARS-COV-2 nucleocapsid protein. The affinity of this antibody binds to SARS-COV-2 nucleocapsid protein is about 2.7 nM.

## Formulation

PBS(pH7.2)

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Human IgG


## Applications

ELISA

## Dilutions

ELISA

## Validations



SDS-PAGE of Anti-N protein Human mAb