

Anti-Nprotein Human mAb



Catalog Number: 166658

Product name

Anti-Nprotein Human mAb

Specificity

Human

Antibody description

human monoclonal antibody (IgG) against SARS-COV-2 nucleocapsid protein. The affinity of this antibody binds to SARS-COV-2 nucleocapsid protein is about 0.66 nM.

Formulation

PBS(pH7.2)

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Human IgG


Applications

ELISA

Dilutions

ELISA

Validations



SDAS-PAGE of Anti-N protein Human mAb