

# Anti-Phospho-EGFR (Tyr1068) Rabbit antibody



Catalog Number: 166670

## Product name

Anti-Phospho-EGFR (Tyr1068) Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human

## Applications

WB

## Antibody description

Phospho-EGFR (Tyr1068) Rabbit polyclonal antibody

## Dilutions

WB: 1/500

## Preparation

Antigen: CDSFLQRP-YSSDPTGA

## Formulation

10mM PBS with 0.02% sodium azide 50% glycerol.


## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Validations



Western blot analysis of extracts of A549 cells, untreated or EGF-stimulated, using Phospho-EGF Receptor (Tyr1068) Rabbit pAb (upper, 1:500 diluted) and EGF Receptor Antibody 201012-6H11 (lower, 1:500 diluted).