

Anti-p38 Rabbit antibody



Catalog Number: 166713

Product name

Anti-p38 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

p38 Rabbit polyclonal antibody

Preparation

Antigen: The antiserum was produced against synthesized peptide derived from human p38 MAPK. AA range:147-196

Formulation

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

IgG

Applications

ICC/IF, WB, IHC-P, ELISA

Dilutions


IF: 1/50-200

WB: 1/500 - 1/2000

IHC: 1/100 - 1/300

ELISA: 1/10000

Validations



Western blot detection of p38 in T47D,3T3,COS7,C6 and HeLa cell lysates using p38 Rabbit pAb (1:5000 diluted).Predicted band size:41KDa.Observed band size:41KDa.