

# Anti-BRCC36 Rabbit antibody



Catalog Number: 166736

## Product name

Rabbit IgG

Anti-BRCC36 Rabbit antibody

## Applications

## Specificity

WB, FC, IP

Human, Mouse, Rat

## Dilutions

## Antibody description

WB: 1/2000-1/10000

BRCC36 Rabbit polyclonal antibody

FC: 1/20-1/100

## Preparation

IP: 1/20-1/50

Antigen: A synthetic peptide of human BRCC36

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage


Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

## Validations



Western blot detection of BRCC36 in Rat Brain, C6, 3T3, HeLa cell lysates using BRCC36 Rabbit pAb(1:1000 diluted). Predicted band size: 36kDa. Observed band size: 36kDa.