

Anti-Annexin VII Rabbit antibody



Catalog Number: 166755

Product name

Anti-Annexin VII Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Annexin VII Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human Annexin VII

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, FC


Dilutions

WB: 1/2000


IHC: 1/50

FC: 1/20

Validations



Immunohistochemistry of Annexin VII in paraffin-embedded Human tonsil using Annexin VII Rabbit



Western blot detection of Annexin VII in K562,Rat Brain cell lysates using Annexin VII Rabbit pAb(1:1000 diluted).Predicted band size:53kDa.Observed band size:53kDa.