

Anti-COX2/Cyclooxygenase 2 Rabbit antibody



Catalog Number: 166838

Product name

Rabbit IgG

Anti-COX2/Cyclooxygenase 2 Rabbit antibody

Applications

Specificity

WB, ICC/IF, FC, IP

Mouse

Dilutions

Antibody description

WB: 1/1000

COX2/Cyclooxygenase 2 Rabbit polyclonal antibody

ICC/IF: 1/20

Preparation

FC: 1/100

Antigen: Recombinant protein of mouse COX2

IP: 1/20

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type



Western blot detection of COX2 in 3t3 cell lysates using COX2 Rabbit pAb(1:1000 diluted).Predicted band size:69KDa.Observed band size:74KDa.