

Anti-Stromal interaction molecule 1 Rabbit antibody



Catalog Number: 166858

Product name

Anti-Stromal interaction molecule 1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Stromal interaction molecule 1 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Stromal interaction molecule 1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000-1/5000

Validations



Western blot detection of Stromal interaction molecule 1 in MCF7,C6,3T3,HeLa cell lysates using Stromal interaction molecule 1 Rabbit pAb(1:1000 diluted).Predicted band size:77KDa.Observed band size:77KDa.