

Anti-p21 Rabbit antibody



Catalog Number: 166862

Product name

Anti-p21 Rabbit antibody

Specificity

Human, Mouse

Antibody description

p21 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human p21

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of p21 in MCF-7 cell lysates using p21 Rabbit pAb(1:1000 diluted).Predicted band size:18KDa.Observed band size:21KDa.