

Anti-p16 ARC Rabbit antibody



Catalog Number: 166891

Product name

Anti-p16 ARC Rabbit antibody

WB: 1/500-1/2000

IHC: 1/50-1/100

Specificity

Human, Mouse, Rat

ICC/IF: 1/50-1/100

IP: 1/50

FC: 1/50

Antibody description

p16 ARC Rabbit polyclonal antibody

Validations

Preparation

Antigen: A synthesized peptide derived from human p16 ARC

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal


Ig Type

Rabbit IgG


Applications

WB, IHC-P, ICC/IF, IP, FC

Dilutions



Immunohistochemical analysis of p16 ARC in paraffin-embedded human colon tissue using p16 ARC antibody at dilution 1/50.



ARC Antibody at dilution 1/50.

Western blot detection of p16 ARC in HepG2,293T,A549,HeLa cell lysates using p16 ARC Rabbit pAb(1:1000 diluted).Predicted band size:16kDa.Observed band size:16kDa.