

Anti-CD147 Rabbit antibody



Catalog Number: 166895

Product name

Anti-CD147 Rabbit antibody

Specificity

Mouse

Antibody description

CD147 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of mouse CD147

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/1000


IHC: 1/20-1/500

ICC/IF: 1/50

FC: 1/50

IP: 1/20

Validations



Western blot detection of CD147/ Basigin in HeLa,C2C12 cell lysates using CD147/ Basigin Rabbit pAb(1:1000 diluted).Predicted band size:42KDa.Observed band size:38-58KDa.