

# Anti-Calnexin Rabbit antibody



Catalog Number: 166910

## Product name

Anti-Calnexin Rabbit antibody

## Specificity

Human, Rat

## Antibody description

Calnexin Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human Calnexin

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, FC

## Dilutions

WB: 1/1000

FC: 1/50

## Validations



Western blot detection of Calnexin in K562 cell lysates using Calnexin Rabbit pAb(1:1000 diluted). Predicted band size:68kDa.Observed band size:90kDa.