

# Anti-SH2D1A/SAP Rabbit antibody



Catalog Number: 166916

## Product name

Anti-SH2D1A/SAP Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

SH2D1A/SAP Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human SH2D1A/SAP

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/1000-1/2000

## Validations



Western blot detection of SH2D1A/SAP in MCF7,C6 cell lysates using SH2D1A/SAP Rabbit pAb(1:1000 diluted).Predicted band size:14KDa.Observed band size:14KDa.