

# Anti-Acetyl-p53 (Lys370) Rabbit antibody



Catalog Number: 166936

## Product name

Anti-Acetyl-p53 (Lys370) Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Acetyl-p53 (Lys370) Rabbit polyclonal antibody

## Preparation

Antigen: A synthesized peptide derived from human p53 (acetyl K370)

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, IP


## Dilutions

WB: 1/500-1/2000

ICC/IF: 1/50-1/200

IP: 1/50

## Validations



Western blot analysis of p53 (acetyl K370) expression in HeLa cell lysate treated with Trichostatin A.