

Anti-PP2A alpha/beta Rabbit antibody



Catalog Number: 166945

Product name

Anti-PP2A alpha/beta Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

PP2A alpha/beta Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human PP2A alpha

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/1000

ICC/IF: 1/20-1/100

FC: 1/20-1/100

IP: 1/50

Validations



Western blot detection of PP2A alpha/beta in K562, Rat Brain, C6, 3T3 cell lysates using PP2A alpha/beta Rabbit pAb(1:1000 diluted). Predicted band size: 36kDa. Observed band size: 36kDa.