

# Anti-COX2/Cyclooxygenase 2 Rabbit antibody



Catalog Number: 166950

## Product name

Rabbit IgG

Anti-COX2/Cyclooxygenase 2 Rabbit antibody

## Applications

## Specificity

WB, ICC/IF, FC, IP

Human, Mouse

## Dilutions

## Antibody description

WB: 1/1000

COX2/Cyclooxygenase 2 Rabbit polyclonal antibody

ICC/IF: 1/50

FC: 1/100

## Preparation

IP: 1/20

Antigen: Recombinant protein of human COX2/Cyclooxygenase 2

## Validations

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type



Western blot detection of COX2/Cyclooxygenase 2 in 3T3 cell lysates using COX2/Cyclooxygenase 2 Rabbit pAb(1:1000 diluted). Predicted band size:69KDa.Observed band size:74KDa.